



Recombinant Mouse ATP synthase F(1) complex catalytic subunit beta, mitochondrial (Atp5f1b)

Product Code	CSB-MP002350MO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P56480
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	AAQASAAPKAGTATGRIVAVIGAVVDVQFDEGLPPILNALEVQGRDSRLVLEVA QHLGESTVRTIAMDGTEGLVRGQKVLDSGAPIKIPVGPETLGRIMNVIGEPIDE RGPIKTKQFAPIHAEAPEFIEMSVEQEILVTGIKVVDLLAPYAKGGKIGLFGGAG VGKTVLIMELINNVAKAHGGYSVFAGVGERTREGNDLYHEMIESGVINLKDATS KVALVYGQMNEPPGARARVALTGLTVAEYFRDQEGQDVLLFIDNIFRFTQAGS EVSALLGRIPSAVGYPQLATDMGTMQERITTTKKGSITSVQAIYVPADDLTDP APATTFAHLDATTVLSRAIAELGIYPAVDPLDSTSRIMDPNIVGNEHYDVARGV QKILQDYKSLQDIILGMDLSEEDKLTVSRARKIQRFLSQPFQVAEVFTGHMG KLVPLKETIKGFQQILAGEYDHLPEQAFYVMVGPIEEAVAKADKLAEHGS
Research Area	Others
Source	Mammalian cell
Target Names	Atp5f1b
Protein Names	Recommended name: ATP synthase subunit beta, mitochondrial EC= 3.6.3.14
Expression Region	47-529aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.