



Recombinant Human Pterin-4-alpha-carbinolamine dehydratase (PCBD1)

Product Code	CSB-MP017514HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P61457
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAGKAHRLSAEERDQLLPNLRAVGVWNELEGRDAIFKQFHFKDFNRAFGFMTR VALQAEKLDHHPWFNVYKVVHITLSTHECAGLSERDINLASFIEQVAVSMT
Research Area	Metabolism
Source	Mammalian cell
Target Names	PCBD1
Protein Names	Recommended name: Pterin-4-alpha-carbinolamine dehydratase Short name= PHS EC= 4.2.1.96 Alternative name(s): 4-alpha-hydroxy-tetrahydropterin dehydratase Dimerization cofactor of hepatocyte nuclear factor 1-alpha Short name=
Expression Region	2-104aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes pterin-4 alpha-carbinolamine dehydratase, an enzyme involved in phenylalanine hydroxylation. A deficiency of this enzyme leads to hyperphenylalaninemia. The enzyme regulates the homodimerization of the transcription factor hepatocyte nuclear factor 1 (HNF1).
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.