




Recombinant Human C-X-C motif chemokine 10 protein (CXCL10)

Product Code	CSB-EP006240HU
Relevance	Chotactic for monocytes and T-lymphocytes. Binds to CXCR3.
Abbreviation	Recombinant Human CXCL10 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P02778
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	VPLSRTVRCTCISISNQPVNPRSLKLEIIPASQFCPRVEIATMKKKGEKRCLNP ESKAIKNLLKAVSKERSKRSP
Research Area	Immunology
Source	E.coli
Target Names	CXCL10
Expression Region	22-98aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	15.5 kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.