



Recombinant Human Estrogen receptor (ESR1), partial

Product Code	CSB-MP007830HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03372
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	SGMALLHQIQGNELEPLNRPQLKIPLERPLGEVYLDSSKPAVYNYPEGAAAYEF NAAAAANAQVYGQTGLPYGPGSEAAAFGSNGLGGFPPLNSVSPSPLMLLHPP PQLSPFLQPHGQQVPYYLENESGYTVREAGPPAFYRPNSDNRRQGGRRERL ASTNDKGSMAKESAKETRYCAVCNDYASGYHYGVWSCEGCKAFFKRSIQGH NDYMCPATNQCTIDKNRRKSCQACRLRKCYEVGMMKGGIRKDRRGGRMLKH KRQRDDGEGRGEVGSAGDMRAANLWPSPLMIKRSKKNLALSALTADQMVSA LLDAEPPILYSEYDPTTRPFSEASMMGLLTNLADRELVHMINWAKRVPGFVDLTL HDQVHLLLECAWLEILMIGLVWRSMHPGKLLFAPNLLLDRNQGKCVVEGMVEIF DMLLATSSRFRMMNLQGEFVCLKSIILLNSGVYTFLSSTLKSLEEKDHIHRVLD KITDTLIHLMAKAGLTLQQQHQRRLAQLLLILSHIRHMSNKGMEHLYSMKCKNVV PLYDLLLEMLDAHRLHAPTSRGGASVEETDQSHLATAGSTSSHSLQKYYITGE AEGFPATV
Research Area	Transcription
Source	Mammalian cell
Target Names	ESR1
Protein Names	Recommended name: Estrogen receptor Short name= ERA Alternative name(s): ER-alpha Estradiol receptor Nuclear receptor subfamily 3 group A member 1
Expression Region	10-595aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Target Details	This gene encodes an estrogen receptor, a ligand-activated transcription factor composed of several domains important for hormone binding, DNA binding, and activation of transcription. The protein localizes to the nucleus where it may form a homodimer or a heterodimer with estrogen receptor 2. Estrogen and its receptors are essential for sexual development and reproductive function, but also play a role in other tissues such as bone. Estrogen receptors are also



involved in pathological processes including breast cancer, endometrial cancer, and osteoporosis. Alternative splicing results in several transcript variants, which differ in their 5 UTRs and use different promoters.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.