




Recombinant Human Protein Wnt-2 (WNT2)

Product Code	CSB-BP026133HU
Abbreviation	Recombinant Human WNT2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P09544
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	SWWYMRATGGSSRVMCDNVPGLVSSQRQLCHRHPDVMRAISQGVAEWTAE CQHQRQHRWNCNTLDRDHSLFGRVLLRSSRESAFVYAISSAGVVFAITRACS QGEVKSCSCDPKKMGSAKDSKGIFDWGGCSDNIDYGIFARAFVDAKERKKGK DARALMNLHNNRAGRKAVKRFLKQECKCHGVSGSCTLRTCWLAMADFRKTG DYLRKYNGAIQVVMNQDGTGFTVANERFKKPTKNDLVYFENSPDYCIRDRE AGSLGTAGRVCNLTSRGMDSCCEVMCCGRGYDTSHVTRMTKCGCKFHWCCA VRCQDCLEALDVHTCKAPKNADWTTAT
Research Area	Cancer
Source	Baculovirus
Target Names	WNT2
Protein Names	Recommended name: Protein Wnt-2 Alternative name(s): Int-1-like protein 1 Int-1-related protein Short name= IRP
Expression Region	26-360aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	41.5
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.