



Recombinant Mouse Glucose-6-phosphate isomerase (Gpi)

Product Code	CSB-MP009717MO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P06745
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	AALTRNPQFQKLEWHRANSANLKLRELFADPERFNNFSLNLTNHHGILVD YSKNLVNKEVMQMLVELAKSRGVEAARDNMFSGSKINYTENRAVLHVALRNR SNTPIKVDGKDVMPENVNRVLDKMKSFQVRSGDWKGYTGKSITDIINIGIGGS DLGPLMVTEALPKPYSKGGPRVWFVSNIDGTHIAKTLASLSPETSLFIIASKTFTT QETITNAETAKEWFLEAAKDPSAVAKHFVALSTNTAKVKEFGIDPQNMFEFWD WVGGRYSLWSAIGLSIALHVGFDFHFEQLLSGAHWMDQHFLKTPLEKNAPVLLA LLGIWYINCYGCETHALLPYDQYMHRFAAYFQQGDMESNGKYITKSGARVDH QTGPIVWGEPTNGQHAFYQLIHQGTKMIPCDFLIPVQTQHPHPIRGLHHKILLA NFLAQTEALMKGKLPPEARKELQAAGKSPEDLEKLLPHKVFEGNRPTNSIVFT KLTPFILGALIAMYEHKIFVQGIMWDINSFDQWGVELGKQLAKKIEPELEGSSAV TSHDSSTNGLISFIKQQRDTKLE
Research Area	Others
Source	Mammalian cell
Target Names	Gpi
Protein Names	Recommended name: Glucose-6-phosphate isomerase Short name= GPI EC= 5.3.1.9Alternative name(s): Autocrine motility factor Short name= AMF Neuroleukin Short name= NLK Phosphoglucose isomerase Short name=
Expression Region	2-558aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene belongs to the GPI family whose members encode multifunctional phosphoglucose isomerase proteins involved in energy pathways. This protein is a dimeric enzyme that catalyzes the reversible isomerization of glucose-6-phosphate and fructose-6-phosphate. The protein functions in different capacities inside and outside the cell. In the cytoplasm, the gene product is involved in glycolysis and gluconeogenesis, while outside the cell it functions as



a neurotrophic factor for spinal and sensory neurons. Defects in this gene are the cause of nonspherocytic hemolytic anemia and a severe enzyme deficiency can be associated with hydrops fetalis, immediate neonatal death and neurological impairment.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.