




Recombinant Human Leukocyte immunoglobulin-like receptor subfamily A member 5 (LILRA5), partial

Product Code	CSB-MP012937HU
Relevance	May play a role in triggering innate immune responses. Does not seem to play a role for any class I MHC antigen recognition.
Abbreviation	Recombinant Human LILRA5 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A6NI73
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	GNLSKATLWAEPGSVISRGN SVTIRCQGTLEAQEYRLVKEGSPEPWDTQNPL EPKNKARFSIPSMTEHHAGRYRCYYYSPAGWSEPSDPLELVVTGFY NKPTLS ALPSPVVTSGENVTLQCGSRLRFDRFILTEEGDHKLSWTLDSQLTPSGQFQAL FPVGPVTPSHRWMLRCYGSRRHILQVWSEPSDLLEIPVSGAADNLSPSQNKS DSGTASHLQDYAVENLIR
Research Area	Immunology
Source	Mammalian cell
Target Names	LILRA5
Expression Region	42-268aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-Myc-tagged
Mol. Weight	29.2kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.