



Recombinant Human Angiopoietin-2 (ANGPT2), partial

Product Code	CSB-BP001707HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O15123
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	YNNFRKSMDSIGKKQYVQVQHGSCSYTFLLPEMDNCRSSSSPYVSNVQRDA PLEYDDSVQRLQVLENIMENNTQWLMKLENIQDNMCKEMVEIQQNAVQNQT AVMIEIGTNLLNQTAEQTRKLT DVEAQVLNQTTTRLELQLEHSLSTNKLEKQILD QTSEINKLQDKNSFLEKKVLAMEDKHIIQLQSIKEEKDQLQVLVSKQNSIIEELEK KIVTATVNNVSLQKQQHDLMETVNNLLTMMSTSNSAKDPTVAKEEQISFRDCA EVFKSGHTTNGIYTLTFPNSTEEIKAYCDMEAGGGGWTIIQRREDGSVDFQRT WKEYKVGFGNPSGEYWLGNEFVSQLTNQQRYVLKIHLDWEGNEAYSLEYH FYLSSEELNYRIHLKGLTGTAGKISSISQPGNDFSTKGDNDKCICKCSQMLTG GWWFDACGPSNLNGMYYPQRQNTNKFNGIKWYYWKGSGYSLK
Research Area	Cardiovascular
Source	Baculovirus
Target Names	ANGPT2
Protein Names	Recommended name: Angiopoietin-2 Short name= ANG-2
Expression Region	19-485aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Target Details	This protein is an antagonist of angiopoietin 1 (ANGPT1) and endothelial TEK tyrosine kinase (TIE-2, TEK). The encoded protein disrupts the vascular remodeling ability of ANGPT1 and may induce endothelial cell apoptosis. Three transcript variants encoding three different isoforms have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.