



Recombinant Rat Angiopoietin-1 (Angpt1)

Product Code	CSB-BP001706RA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O35460
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	SNQRRSPENGGRRYRNRIQHGQCAYTFILPEHDGNCRESATEQYNTNALQRDA PHVETDFSSQKLQHLEHVMENYTQWLQKLENYIVENMKSEMAQIQQNAVQNH TATMLEIGTSLLSQTAEQTRKLTVDVETQVLNQTSLRLEIQLLNSLSTYELEKQLL QQTNEILKIQEKNSLLEHKILEMEGKHKEELDTLKEEKENLQGLVTRQTFIIQELE KQLSRATSNNSVLQKQQLLELMDTVHNLVSLCTKEVLLKGGKREEEKPF RDCA DVYQAGFNKSGIYTIYFNNMPEPKKVF CNMDVNEGGWTVIQHREDGSLDFQR GWKEYKMGFGNPSGEYWLGNEFIFAITSQRQYMLRIELMDWEGNRAYSQYD RFHIGNQKQNYRLYLKGHTGTAGKQSSLILHGADFSTKDADNDNCMCKCALM LTGGWWFDACGPSNLNGMFYTAGQNHGKLNKWHYFKGPSYSLRSTTMMI RPLDF
Research Area	Others
Source	Baculovirus
Target Names	Angpt1
Protein Names	Recommended name: Angiopoietin-1 Short name= ANG-1
Expression Region	20-497aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	Angiopoietins are proteins with important roles in vascular development and angiogenesis. All angiopoietins bind with similar affinity to an endothelial cell-specific tyrosine-protein kinase receptor. This protein is a secreted glycoprotein that activates the receptor by inducing its tyrosine phosphorylation. It plays a critical role in mediating reciprocal interactions between the endothelium and surrounding matrix and mesenchyme. The protein also contributes to blood vessel maturation and stability, and may be involved in early development of the heart.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.