



Recombinant Mouse Paired immunoglobulin-like receptor B (Pirb), partial

| | |
|--------------------------|---|
| Product Code | CSB-MP012941MO |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P97484 |
| Storage Buffer | Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0 |
| Product Type | Recombinant Proteins |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | Greater than 85% as determined by SDS-PAGE. |
| Sequence | RRRHRGKFRKDVQKEKDLQLSSGAEEPITRKGELQKRPNPAAATQEESLYAS VEDMQTEDGVELNSWTPPEEDPQGETYAQVKPSRLRKAGHVSPSVMSREQL NTEYEQAEEGQGANNQAAESGESQDVTYAQLCSRTLRLQGAAASPLSQAGEA PEEPSVYATLAAARPEAVPKDMEQ |
| Research Area | Others |
| Source | Mammalian cell |
| Target Names | Pirb |
| Protein Names | Recommended name: Leukocyte immunoglobulin-like receptor subfamily B member 3 Short name= LIR-3 Short name= Leukocyte immunoglobulin-like receptor 3Alternative name(s): Cell-surface glycoprotein p91 Paired immunoglobulin-like re |
| Expression Region | 664-841aa |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Partial |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |