



Recombinant Helicobacter pylori Urease subunit alpha (ureA)

Product Code	CSB-BP320452HUV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P14916
Product Type	Recombinant Protein
Immunogen Species	Helicobacter pylori (strain ATCC 700392 / 26695) (Campylobacter pylori)
Purity	≥85% (SDS-PAGE)
Sequence	MKLTPEKELDKLMLHYAGELAKKRKEKGIKLNYVEAVALISAHIMEEARAGKKTAAELMQEGRTLLKPDDVMDGVASMIHEVGIAMFPDGTKLVTVHTPIEANGKLVPGELFLKNEDITINEGKKAHSVVKVKNVGDPRVQIGSHFFHFEVNRCLDFDREKTFGKRLDIASGTAVRFEPGEEKSVELIDIGGNRRIFGFNALVDRQADNESKKIALHRAKERGFHGAKSSDNYVKTIKE
Research Area	Signal Transduction
Source	Baculovirus
Target Names	ureA
Protein Names	Recommended name: Urease subunit alpha EC= 3.5.1.5 Alternative name(s): Urea amidohydrolase subunit alpha
Expression Region	1-238aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.