



Recombinant *Vibrio cholerae* serotype O1 Zona occludens toxin (zot)

Product Code	CSB-EP334425VEX-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38442
Product Type	Recombinant Protein
Immunogen Species	<i>Vibrio cholerae</i> serotype O1 (strain ATCC 39315 / El Tor Inaba N16961)
Purity	>85% (SDS-PAGE)
Sequence	MSIFIHGAPGSYKTS GALWLRLLPAIKSGRHIITNVRGLNLERMAKYLKMDVS DISIEFIDTDHPDGRLT MARFWHWARKDAFLFIDECGRIWPPRLT VTNLKALDT PPDLVAEDRPESFEVAFDMHRHHGWDICL TTPNIAKVHNMIREAAEIGYRHFN RATVGLGAKFTLTTHDAANS GQMDSHALTRQVKKIPSPIFKMYASTTTGKARD TMAGTALWKDRKILFLFGMVFLMFSYSFYGLHDNPIFTGGNDATIESEQSEPO SKATVGNVAVGSKAVAPASFGFCIGRLCVQDGFVTVGDERYRLVDNLDIPYRGL WATGHHIYKDTLTVFFETESGSVPT ELFASSYRYKVLPLPDFNHFVVFDTFAAQ ALWVEVKRGLPIKTENDKKGLNSIF
Research Area	Others
Source	E.coli
Target Names	zot
Protein Names	Recommended name: Zona occludens toxin Alternative name(s): Zonular occludens toxin Short name= Zot
Expression Region	1-399aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.