



Recombinant *Xenopus laevis* Protein Wnt-8 (wnt8)

Product Code	CSB-MP340936XBE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P28026
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	AWSVNNFL MTGPKAYLTY SASVAVGAQN GIEECKYQFA WERWNCPEST LQLATHNGLR SATRETSFVH AISSAGVMYT LTRNCSMGDF DNCGCDDSRN GRIGGRGWVW GGCSDNAEFG ERISKLFVDG LETGQDARAL MNLHNNEAGR LAVKETMKRT CKCHGISGSC SIQTCWLQLA EFRDIGNHLK IKHDQALKLE MDKRKMRSNGN SADNRGAIAD AFSSVAGSEL IFLEDSPDYC LKNISLGLQG TEGRECLQSG KNLSQWERRS CKRLCTDCGL RVEEKKTEII SSCNCKFHWC CTVKCEQCKQ VVIKHFCARR ERDSNMLNTK RKNRGHRR
Source	Mammalian cell
Target Names	wnt8
Protein Names	Recommended name: Protein Wnt-8 Short name= XWnt-8
Expression Region	23-358
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.