



Recombinant Danio rerio D-2-hydroxyglutarate dehydrogenase, mitochondrial (d2hgdh)

Product Code	CSB-BP374600DIL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A1L258
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	ANPAN PSAAPPRLPF SRVTQEDLSF FRALLPGRTI TDPDLLKSSN VDWLKTVQGS SDVLLRPKTT EGVSQILRYC NERNLAVCPQ GGNTGLVGGG VPVFDEIILS TSLMNQVFAF DNISGILTCQ AGCVLENLSH YLEERDFIMP LDLGAKGSCH IGGNVSTNAG GLRLLRYGSL RGTVLGLEVV LADGHVLNCL ATLRKDNTGY DLKQLFIGSE GTLGVITAVS ILCPRKPKAV NVAFLGCSSF QQLLETQCC RGMLGEILSA FEFLDASCMN LLEKHLKLTN PITECPFYIV IETAGSNATH DEEKLHQFLE EVMTSSLVTD GTVATEATKI KALWSLRERV TEALTHEGYT YKYDISLPVE KIYDLVQDMR RHLGGMAKNV VGYGHVGDGN LHLNITSPSK DFDLLAAIEP YVYEWTSQWK GSISAEHGLG LKKRNYIYYS KPSEAVALMG SIKAMLDPKG ILNPKYKTLPD NIN
Source	Baculovirus
Target Names	d2hgdh
Protein Names	Recommended name: D-2-hydroxyglutarate dehydrogenase, mitochondrial EC=1.1.99.-
Expression Region	56-533
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.