



Recombinant Callithrix jacchus Interleukin-10 (IL10)

Product Code	CSB-MP611872CYL
Abbreviation	IL10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0Z972
Product Type	Recombinant Protein
Immunogen Species	Callithrix jacchus (White-tufted-ear marmoset)
Purity	≥85% (SDS-PAGE)
Sequence	SPGQGTQSENSCTHFPGLPHMLRELRVAFSRVKTFQKDKQLDSMLLKESL LEDFKGYLGQCQALSEMIQFYLEEVMQAENHDPDIKEHVNSLGEKLTFRRLR RRCHRFLPCENKSKAVVQVKNVSKLQEKGIYKAMSEFDIFIDYIEAYMTMKAQ N
Source	Mammalian cell
Target Names	IL10
Protein Names	Recommended name: Interleukin-10 Short name= IL-10 Alternative name(s): Cytokine synthesis inhibitory factor Short name= CSIF
Expression Region	19-178aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.