



Recombinant Human Dihydropyrimidinase-related protein 2 (DPYSL2)

Product Code	CSB-BP614975HU
Abbreviation	DPYSL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q16555
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSYQGKKNIPRITSDRLLIKGGKIVNDDQSFYADIYMEDGLIKQIGENLIVPGGV KTIEAHSRMVIPGGIDVHTRFQMPDQGMTSADDFQGTKAALAGGTTMIIDHV VPEPGTSLAAFDQWREWADSKSCCDYSLHVDISEWHKGIQEEEMEALVKDHG VNSFLVYMAFKDRFQLTDCQIYEVLSVIRDIGAIQVHAENGDIIEEQQRILDLG ITGPEGHVLSRPEEVEAEAVNRAITIANQTNCPYITKVMSSAEVIAQARKKG TVVYGEPIASLGTGSHYWSKNWAKAAAFVTSPPLSPDPTTPDFLNSLLSCG DLQVTGSAHCTFNQAQKAVGKDNFTLIPEGTNGTEERMSVIWDKAVVTGKMD ENQFVAVTSTNAAKVFNLYPRKGRIAVGSDADLVIWDPDSVKTISAKTHNSSLE YNIFEGMECRGSPLVVISQGKIVLEDGTLHVTEGSGRYIPRKPFDFVYKRIKA RSRLAELRGVPRGLYDGPVCEVSVTPKTVTPASSAKTSPAKQQAPPVRNLHQ SGFSLSGAQIDDNIPRRTTQRIVAPPGGRANITSLG
Research Area	Developmental Biology
Source	Baculovirus
Target Names	DPYSL2
Protein Names	Recommended name: Dihydropyrimidinase-related protein 2 Short name= DRP-2 Alternative name(s): Collapsin response mediator protein 2 Short name= CRMP-2 N2A3 Unc-33-like phosphoprotein 2 Short name= ULIP-2
Expression Region	1-572aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.