



# Recombinant Human Inactive tyrosine-protein kinase 7 (PTK7), partial

<b>Product Code</b>	CSB-MP622651HU
<b>Abbreviation</b>	PTK7
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q13308
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	AIVFIKQPSSQDALQGRRALLRCEVEAPGPVHVYWLLDGAPVQDTERRFAQG SSLSFAAVDRLQDSGTFQCVARDDVTGEEARSANASFNKWKIEAGPVVLKHPA SEAEIQPQTQVTLRCHIDGHRPTYQWFRDGTPLSDGQSNHTVSSKERNLTL RPAGPEHSGLYSCCAHSAFGQACSSQNFTLSIADESFARVVLAPQDVVVARY EEAMFHCQFSAQPPPSLQWLFEDETPITNRSRPPHLRRATVFANGSLLLTQVR PRNAGIYRCIGQGQRGPPIILEATLHLAEIEDMPLFEPVFTAGSEERTVCLPPK GLPEPSVWWEHAGVRLPTHGRVYQKGHELVLANIAESDAGVYTCHAANLAG QRRQDVNITVATVPSWLKKPQDSQLEEGKPGYLDCLTQATPKPTVVWYRNQ MLISEDSRFEVFKNGTLRINSVEVYDGTWYRCMSSTPAGSIEAQRVQVLEKL KFTPPPQPQCMEFDKEATVPCSATGREKPTIKWERADGSSLPEWVTDNAGT LHFARVTRDDAGNYTCIASNGPQGQIRAHVQLTVAVFITFKVEPERTTVYQGH TALLQCEAQGDPKPLIQWKGKDRILDPTKLGPRMHIFQNGSLVIHDVAPEDSG RYTCIAGNSCNIKHTEAPLYVVDKPVPEESEGGSPPPYKMIQT
<b>Research Area</b>	Cell Adhesion
<b>Source</b>	Mammalian cell
<b>Target Names</b>	PTK7
<b>Protein Names</b>	Recommended name: Inactive tyrosine-protein kinase 7 Alternative name(s): Colon carcinoma kinase 4 Short name= CCK-4 Protein-tyrosine kinase 7 Pseudo tyrosine kinase receptor 7 Tyrosine-protein kinase-like 7
<b>Expression Region</b>	31-704aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.