



# Recombinant Human Inactive tyrosine-protein kinase 7 (PTK7), partial

<b>Product Code</b>	CSB-MP622651HU
<b>Abbreviation</b>	PTK7
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q13308
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 85% as determined by SDS-PAGE.
<b>Sequence</b>	AIVFIKQPSSQDALQGRRALLRCEVEAPGPVHVYWLLDGAPVQDTERRFAGG SSLSFAAVDRLQDSGTFQCVARDDVTGEEARSANASFNIKWIEAGPVVLKHPA SEAEIQPQTQVTLRCHIDGHPRPTYQWFRDGTPLSDGQSNHTVSSKERNLTL RPAGPEHSGLYSCCAHSAFGQACSSQNFTLSIADES FARVVLAPQDVVVARY EEAMFHCQFSAQPPSLQWLFEDETPITNRSRPPHLRRATVFANGSLLLTQVR PRNAGIYRCIGQGQRGPPPIILEATLHLAEIEDMPLFEP RVFTAGSEERTVCLPPK GLPEPSVWWEHAGVRLPTHGRVYQKGHELVLANIAESDAGVYTCHAANLAG QRRQDVNITVATVPSWLKPKQDSQLEEGKPGYLDCLTQATPKPTVVWYRNQ MLISED SRFEVFKNGTLRINSVEVYDGTWYRCMSSTPAGSIEAQRVQVLEKL KFTPPPQPQQCMEFDKEATVPCSATGREKPTIKWERADGSSSLPEWVTDNAGT LHFARVTRDDAGNYTCIASNGPQGQIRAHVQLTVAVFITFKVEPERTTVYQGH TALLQCEAQGDPKPLIQWKGKDRILDPTKLGPRMHIFQNGSLVIHADVAPEDSG RYTCIAGNSCNIKHTEAPLYVVDKPVPEESEGGSPPPYKMIQT
<b>Research Area</b>	Cell Adhesion
<b>Source</b>	Mammalian cell
<b>Target Names</b>	PTK7
<b>Protein Names</b>	Recommended name: Inactive tyrosine-protein kinase 7 Alternative name(s): Colon carcinoma kinase 4 Short name= CCK-4 Protein-tyrosine kinase 7 Pseudo tyrosine kinase receptor 7 Tyrosine-protein kinase-like 7
<b>Expression Region</b>	31-704aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.