



Recombinant Human Protein Wnt-10a (WNT10A)

Product Code	CSB-MP884424HU
Abbreviation	WNT10A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9GZT5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPSA PNDILDRLRP PEPVLNANTV CLTLPGLSRR QMEVCVRHPD VAASAIQGIQ IAIHECQHQF RDQRWNCSSL ETRNKIPYES PIFSRGFRES AFAYAIAAAG VVHAVSNACA LGKLGKACGCD ASRRGDEEAF RRKLHRLQLD ALQRGKGLSH GVPEHPALPT ASPGLQDSWE WGGCSPDMGF GERFSKDFLD SREPHRDIHA RMRLHNNRVG RQAVMENMRR KCKCHGTSGS CQLKTCWQVT PEFRTVGALL RSRFHRATLI RPHNRNGGQL EPGPAGAPSP APGAPGPRRR ASPADLVYFE KSPDFCEREP RLDSAGTVGR LCNKSSAGSD GCGSMCCGRG HNILRQTRSE RCHCRFWCC FVVCCECRIT EWVSVCK
Source	Mammalian cell
Target Names	WNT10A
Protein Names	Recommended name: Protein Wnt-10a
Expression Region	36-417
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene is a member of the WNT gene family. It is strongly expressed in the cell lines of promyelocytic leukemia and Burkitt's lymphoma. In addition, it and another family member, the WNT6 gene, are strongly coexpressed in colorectal cancer cell lines. The gene overexpression may play key roles in carcinogenesis through activation of the WNT-beta-catenin-TCF signaling pathway. This gene and the WNT6 gene are clustered in the chromosome 2q35 region.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.