




Recombinant Human Dystonin (DST), partial

Product Code	CSB-EP821036HU
Relevance	Cytoskeletal linker protein. Acts as an integrator of intermediate filaments, actin and microtubule cytoskeleton networks. Required for anchoring either intermediate filaments to the actin cytoskeleton in neural and muscle cells or keratin-containing intermediate filaments to hemidesmosomes in epithelial cells. The proteins may self-aggregate to form filaments or a two-dimensional mesh. Isoform 3: plays a structural role in the assembly of hemidesmosomes of epithelial cells; anchors keratin-containing intermediate filaments to the inner plaque of hemidesmosomes. Required for the regulation of keratinocyte polarity and motility; mediates integrin ITGB4 regulation of RAC1 activity. Isoform 6: required for bundling actin filaments around the nucleus. By similarity3 Publications Isoform 7: regulates the organization and stability of the microtubule network of sensory neurons to allow axonal transport.
Abbreviation	Recombinant Human DST protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q03001
Alias	230 kDa bullous pemphigoid antigen 230/240 kDa bullous pemphigoid antigen Bullous pemphigoid antigen 1 Short name:BPA Short name:Bullous pemphigoid antigen
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MHSSSYSYRSSDSVFSNTTSTRTSLDSNENLLLHCGPTLINSICISFGSESF DGHRLMLQQIANRVQRDSVICEDKLILAGNALQSDSKRLESGVQFQNEAEIAGYI LECENLLRQHVIDVQILIDGKYYQADQLVQRVAKLRDEIMALRNECSSVYSKGR ILTTEQTKLMISGITQSLNSGFAQTLHPSLTSG
Research Area	Signal Transduction
Source	E.coli
Target Names	DST
Expression Region	1-195aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	39.1kDa
Protein Length	Partial of Isoform 3



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.