



Recombinant pig Trypsin

Product Code	CSB-EP018814PI
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P00761
Product Type	Recombinant Protein
Immunogen Species	Sus scrofa (Pig)
Purity	>85% (SDS-PAGE)
Sequence	IVGGYTCAANSIPYQVSLNSGSHFCGGSLINSQWVVSAAHICYKSRIQVRLGEH NIDVLEGNEQFINAAKIITHPNFNGNTLDNDIMLIKLSPPATLNSRVATVSLPRSC AAAGTECLISGWGNTKSSGSSYPSLLQCLKAPVLSOSSCKSSYPGQITGNMIC VGFLEGGKDSCQGDSSGPPVVCNGQLQGIVSWGYGCAQKNKPGVYTKVCNY VNWIQQTIAAN
Research Area	Others
Source	E.coli
Expression Region	9-231aa
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.