



Recombinant *Macaca mulatta* Interleukin-1 alpha (IL1A)

Product Code	CSB-EP011613MOW-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P48089
Product Type	Recombinant Protein
Immunogen Species	<i>Macaca mulatta</i> (Rhesus macaque)
Purity	≥85% (SDS-PAGE)
Sequence	SAPFSFLS NMTYHFIRII KHEFILNDTL NQTIIRANDQ HLTAAAIHNL DEAVKFDMGGA YTSSKDDTKV PVILRISKTKQ LYVSAQDEDQ PVLLKEMPEI NKTITGSETN FLFFWETHGT KNYFISVAHP NLFIATKHNDN WVCLAKGLPS ITDFQILENQ A
Source	<i>E. coli</i>
Target Names	IL1A
Protein Names	Recommended name: Interleukin-1 alpha Short name= IL-1 alpha Alternative name(s): Hematopoietin-1
Expression Region	113-271
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is a member of the interleukin 1 cytokine family. This cytokine is a pleiotropic cytokine involved in various immune responses, inflammatory processes, and hematopoiesis. This cytokine is produced by monocytes and macrophages as a proprotein, which is proteolytically processed and released in response to cell injury, and thus induces apoptosis. This gene and eight other interleukin 1 family genes form a cytokine gene cluster on chromosome 2. It has been suggested that the polymorphism of these genes is associated with rheumatoid arthritis and Alzheimer's disease.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.