



Recombinant Rat Allograft inflammatory factor 1 (Aif1)

Product Code	CSB-MP001490RA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P55009
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	SQSKDLQGG KAFGLLKAQQ EERLDGINKH FLDDPKYSSD EDLQSKLEAF KTKYMEFDLN GNGDIDIMSL KRMLEKLGVP KTHLELKKLI REVSSGSEET FSYSDFLRMM LGKRSAILRM ILMYEEKNKE HQKPTGPPAK KAISELP
Source	Mammalian cell
Target Names	Aif1
Protein Names	Recommended name: Allograft inflammatory factor 1 Short name= AIF-1 Alternative name(s): Ionized calcium-binding adapter molecule 1 MRF-1 Microglia response factor
Expression Region	2-147
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene is induced by cytokines and interferon. Its protein product is thought to be involved in negative regulation of growth of vascular smooth muscle cells, which contributes to the anti-inflammatory response to vessel wall trauma. Three transcript variants encoding different isoforms have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.