



Recombinant Mouse Protein disulfide-isomerase A2 (Pdia2)

Product Code	CSB-MP017719MO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	D3Z6P0
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	QGEPPGGPSEVLPEEPTGEEVPKEDGILVLNHR T LSLALQEHSALMVEFYAPW CGHCKELAPEYSKAAALLAAESAVVTLAKVDGPAEPELTKEFEVVGYP T LKFF QNGNRTNPEEYAGPKTAEGIAEWLRRRVG P SATHLEDEEGVQALMAKWDMV VIGFFQDLQGKDMATFLALAKDALDMTFGFTDQPQLFEKFG LTKD T VVLFKKF DEGRADFPVDKETGLDLGDL SRFLVIHSMHLVTEFNSQTSPKIFA AKILNHLLLF VNQTLAQHRELLTDFREAA P PFRGQVLFVMVDVAADNSHVLN YFGLKAE EAP TLRLIN VETTKKYAPTGVIAITAASVAAFCQAVLHGEIKHYLLS QEIPPDWDQGP VKTLVSKNFEQVAFDETKNV FVKFYAPWCSHCKEMAPAW EALAEKYKDREDI VIAELDATANELEAFSVLGYPTL KFFPAGPDRKVIDYKSTRDLETFSKFLD SGG HLPKEEPKEPAASAPEAQANSTLGPKEEL
Research Area	Others
Source	Mammalian cell
Target Names	Pdia2
Protein Names	Recommended name: Protein disulfide-isomerase A2 EC= 5.3.4.1 Alternative name(s): PDIp
Expression Region	21-527aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.