



# Recombinant Human Long-chain fatty acid transport protein 2 (SLC27A2), partial

<b>Product Code</b>	CSB-EP021534HU-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O14975
<b>Storage Buffer</b>	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	GATLALRTKFSASQFWDDCRKYNVTVIQYIGELRLYLCNSPQKPNDRDHKVRL ALGNGLRGDVWRQFVKRFGDICIYEFYAATEGNIGFMNYARKVGVAVGRVNYL QKKIITYDLIKYDVEKDEPVRDENGVCVRVPKGEVGLLVCKITQLTPFNGYAGA KAQTEKKKLRDVFVKGDLYFNSGDLLMVDHENFIYFHDRVGDTRFWKGENVA TTEVADTVGLVDFVQEVNVYGVHVPDHEGRIGMASIKMKENHEFDGKKLFQHI ADYLPSYARPRFLRIQDTIEITGTFKHRKMTLVEEGFNPAVIKDALYFLDDTAKM YVPMTEDIYNAISAKTLKL
<b>Research Area</b>	Metabolism
<b>Source</b>	E.coli
<b>Target Names</b>	SLC27A2
<b>Protein Names</b>	Recommended name: Very long-chain acyl-CoA synthetase Short name= VLACS Short name= VLCS EC= 6.2.1.-Alternative name(s): Fatty acid transport protein 2 Short name= FATP-2 Fatty-acid-coenzyme A ligase, very long-chain 1
<b>Expression Region</b>	283-620aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.