



Recombinant *Acinetobacter baumannii* Elongation factor Tu (tuf1)

Product Code	CSB-EP025342AUD-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A3M1F6
Product Type	Recombinant Protein
Immunogen Species	<i>Acinetobacter baumannii</i> (strain ATCC 17978 / CIP 53.77 / LMG 1025 / NCDC KC755 / 5377)
Purity	≥85% (SDS-PAGE)
Sequence	MAKAKFERNKPHVNVGTIGHVDHGKTTLTAAIATICAITYGGGEAKDYSQIDSAP EEKARGITINTSHVEYDSPTRHYAHVDCPGHADYVKNMITGAAQMDGAILVCA ATDGPMPQTREHILLSRQVGVPYIIVFLNKCDLVDDEELLELVEMEVRELLSTY DFPGDDTPVIRGSALAALNGEAGPYGEE SVLALVAALDSYIPEPERAIDKAFLM PIEDVFSISGRGTVVTGRVEAGIIVGVEEVEIVGIKDTVKTVTGVEMFRKLLDE GRAGENCGILLRGTREEVQRGQVLAKPGTIKPHTKFDAEVVLSKEEGGRHT PFLNGYRPQFYFRITTDVTGAIQLKEGVEMVMPGDNVEMSVELIHPIAMDPGLR FAIREGGRTVGAGVVAKVTA
Source	E.coli
Target Names	tuf1
Protein Names	Recommended name: Elongation factor Tu Short name= EF-Tu
Expression Region	1-396aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.