



Recombinant *Cryptomeria japonica* Pectate lyase

1

Product Code	CSB-BP322455DYO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P18632
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	<i>Cryptomeria japonica</i> (Japanese cedar) (<i>Cupressus japonica</i>)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	DNPIDSCWRGDSNWAQNRMKLADCAVGFSGSSTMGGKGGDLYTVTNSDDDP VNPAPGTLRYGATRDRPLWIIFSGNMNIKMKMPMYIAGYKTFDGRGAQVYIGN GGPCVFIKRVSNVIIHGLHLYGCSTSVLGNVLINESFGVEPVHPQDGDALTLRT ATNIWIDHNSFSNSSDGLVDVTLSSSTGVTISNNLFFNHHKVMMLLGHDDAYSDD KSMKVTVAFNQFGPNCGQRMPRARYGLVHVANNYDPWTIYAIGGSSNPTIL SEGENSFTAPNESYKKQVTIRIGCKTSSSCSNWVWQSTQDVFYNGAYFVSSGK YEGGNIYTKKEAFNVENGNATPQLTKNAGVLTCSLSKRC
Research Area	Others
Source	Baculovirus
Target Names	Allergen Cry j I
Protein Names	Recommended name: Sugi basic protein Short name= SBPAlternative name(s): Allergen Cry j I Major allergen Cry j 1 Allergen= Cry j 1
Expression Region	22-374aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.