



Recombinant Aspergillus niger 3-phytase A (phyA)

Product Code	CSB-BP333934DOZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34752
Product Type	Recombinant Protein
Immunogen Species	Aspergillus niger
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	ASRNQSS CDTV DQGYQC FSETSHLWGQ YAPFFSLANE SVISPEVPAG CRVTF AQVLS RHGARYPTDS KGKKYSALIE EIQQNATTFD GKYAFLKTYN YSLGADDLTP FGEQELVNSG IKFYQRYESL TRNIVPFIRS SGSSRVIASG KKFIEGFQST KDKPRAQPG QSSPKIDVVI SEASSSNTL DPGTCTVFED SELADTVEAN FTATFVPSIR QRLENDLSGV TLTDEVTYL MDMCSFDTIS TSTVDTKLSP FCDLFTHDEW INYDY LQSLK KYYGHGAGNP LGPTQGVGYA NELIARLTHS PVHDDTSSNH TLDSSPATFP LNSTLYADFS HDNGIISILF ALGLYNGTKP LSTTTVENIT QTDGFSSAWT VPFASRLYVE MMQCQAEQEP LVRVLVNDRV VPLHGCPVDA LGRCTRDSFV RGLSFARSGG DWAECA
Source	Baculovirus
Target Names	phyA
Protein Names	Recommended name: 3-phytase A EC= 3.1.3.8 Alternative name(s): 3 phytase A Myo-inositol hexakisphosphate phosphohydrolase A Myo-inositol-hexaphosphate 3-phosphohydrolase A
Expression Region	24-467
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.