



Recombinant Escherichia coli 1-acyl-sn-glycerol-3-phosphate acyltransferase (plsC)

| | |
|--------------------------|---|
| Product Code | CSB-MP340803ENV |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P26647 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli (strain K12) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MLYIFRLIITVIYSILVCFVFGSIYCLFSPRNPKHVATFGHMFGRLLAPLFGGLKVECR KPTDAESYGNAIYIANHQNNYDMVTASNIVQPPTVTVGKKSLLWIPFFGQLYWL TGNLLIDRNNRTKAHGTTAEVVNHFKRRISIWMFPEGTRSRGRGLLPFKTGAF HAAIAAGVPIIPVVCVSTTSNKINLNRLHNGLVIVEMLPIDVSQYGKDKQVRELAA HCRSIMEQKIAELDKVAEREAAGKV |
| Research Area | Others |
| Source | Mammalian cell |
| Target Names | plsC |
| Protein Names | Recommended name: 1-acyl-sn-glycerol-3-phosphate acyltransferase Short name= 1-AGP acyltransferase Short name= 1-AGPAT EC= 2.3.1.51 EC= 2.3.1.n4 Alternative name(s): Lysophosphatidic acid acyltransferase Short name= |
| Expression Region | 1-245aa |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |