



Recombinant Zaire ebolavirus Pre-small/secreted glycoprotein (GP)

| | |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-MP350772ZAT |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P60171 |
| Product Type | Recombinant Protein |
| Immunogen Species | Zaire ebolavirus (strain Kikwit-95) (ZEBOV) (Zaire Ebola virus) |
| Purity | >85% (SDS-PAGE) |
| Sequence | IPLGVIHN STLQVSDVDK LVCRDKLSST NQLRSVGLNL EGNQVATDVP SATKRWGFRS GVPPKVVNYE AGEWAENCYN LEIKKPDGSE CLPAAPDGIR GFPRCRYVHK VSGTGPCAGD FAFHKEGAFF LYDRLASTVI YRGTTFAEGV VAFLILPQAK KDFSSHPLR EPVNATEDPS SGYYSTTIKY QATGFGTNET EYLFVDNLT YVQLESRFTP QLLQLNETI YTSKRSNTT GKLIWKNPE IDTTIGEWAF WETKKSLEK FAVKSCSQL YQTEPKTSVV RVRRELLPTQ GPTQQLKTTK SWLQKIPLQW FKCTVKEGKL QCRI |
| Source | Mammalian cell |
| Target Names | GP |
| Protein Names | Recommended name: Pre-small/secreted glycoprotein Short name= pre-sGP Cleaved into the following 2 chains: 1. Small/secreted glycoprotein Short name= 2. sGP 3. Delta-peptide |
| Expression Region | 33-364 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |