



Recombinant Mycobacterium tuberculosis Steroid C26-monooxygenase (cyp125)

Product Code	CSB-MP351308MVZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P9WPP1
Product Type	Recombinant Protein
Immunogen Species	Mycobacterium tuberculosis (strain ATCC 25618 / H37Rv)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MSWNHQSVEIAVRRTTVSPNLPPGFDFDTPAIYAERLPVAEFAELRSAAPIW WNGQDPGKGGGFHDGGFWAITKLNDVKEISRHSDFSSYENGVIPRFKNDIA REDIEVQRFVMLNMDAPHHTRLRKIISRGFTPRAVGRLHDELQERAQKIAAEAA AAGSGDFVEQVSCLEPLQAIAGLLGVPQEDRGKLFHWSNEMTGNEDPEYAH DPKASSAELIGYAMKMAEEKAKNPADDIVTQLIQADIDGKLSDDFFGVVML AVAGNETTRNSITQGMMAFAEHPDQWELYKKVRPETADEIVRWATPVTAFAQ RTALRDYELSGVQIKKGQRVVMFYRSANFDEEVFQDPFTFNILRNPNPHVGF GTGAHYCIGANLARMTINLIFNAVADHMPDLKPISAPERLRSGWLNGIKHWQV DYTGRCPPVAH
Source	Mammalian cell
Target Names	cyp125
Protein Names	Recommended name: Steroid C27-monooxygenase EC= 1.14.13.-Alternative name(s): Cholest-4-en-3-one monooxygenase Cytochrome P450 125
Expression Region	1-433aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.