



Recombinant Escherichia phage T7 Single-stranded DNA-binding protein (2.5)

Product Code	CSB-BP366021EEB
Abbreviation	Recombinant Escherichia phage T7 2.5 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03696
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Escherichia phage T7 (Bacteriophage T7)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MAKKIFTSALGTAEPYAYIAKPDYGNEERGFGNPRGVYKVDLTIPNKDPRCQR MVDEIVKCHEEAYAAAVEEYEANPPAVARGKKPLKPYEGDMPFFDNGDGGTTT FKFKCYASFQDKKTKETKHINLVVVDVSKGKKMEDVPIIGGGSKLKVYSLVPYK WNTAVGASVKLQLESVMLVELATFGGGEDDWADEVEENGYVASGSAKASKP RDEESWDEDEESEEAEDEDGDF
Source	Baculovirus
Target Names	2.5
Protein Names	Recommended name: Helix-destabilizing protein Alternative name(s): Single-stranded DNA-binding protein Short name= SSB protein
Expression Region	1-232aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.