



Recombinant Rabies virus Glycoprotein (G), partial

Product Code	CSB-BP389005RIE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A3RM22
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Rabies virus (strain India) (RABV)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	KFPIYTIPDKLGPWSPIDIHHLSCPNNLVEDEGCTNLSGFSYMELKVGYISAIK VNGFT CTGVVTEAETYTNFVG YVTTTTFKRKHFRPTDACRAAYNWKMAGDPRYEESL HNPYPDYH WLRTVKTTKESLVIISPSVADLDPYDKSLHSRVFPSGKCSGITISSTYCSTNHDY TIWMP ENPRLGTSCDIFTNSRGKRASKGGKTCGFVDERGLYKSLKGACKLKL CGVLGL RLMDGTW VAMQTSDETKWCPPDQLVNLHDFRSDEIEHLVVEELVKKREECLDALESIMAT KSVSFRR LSHLRKLVPFGFKAYTIFNKTLMEADAHYKSVRTWNEIIPSKGCLR VGGRCHP HVNGVFF NGIILGPDGHVLIPEMQSSLLQQHMELLESSVIPLMHPLADPSTVFKDGDEAED FVEVHL PDVHKQISGVDLGLPSWGKY
Source	Baculovirus
Target Names	Name:G
Protein Names	Recommended name: Glycoprotein G
Expression Region	20-459
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.