



# Recombinant *Klebsiella pneumoniae* Metallo-beta-lactamase type 2 (blaNDM-1)

<b>Product Code</b>	CSB-BP510513KBG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	C7C422
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Klebsiella pneumoniae</i>
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	GEIRPTIGQQMETGDQRFGDLVFRQLAPNVWQHTSYLDMPGFGAVASNGLIV RDGGRVLVVDTAWTDDQTAQILNWIKQEINLPVALAVVTHAHQDKMGGMDAL HAAGIATYANALSNQLAPQEGMVAAQHSLTFAANGWVEPATAPNFGPLKVFY PGPGHTSDNITVGIDGTDIAFGGGLIKDSKAKSLGNLGDADTEHYAASARAFGA AFPKASMIVMSHSAPDSRAAITHTARMADKLR
<b>Source</b>	Baculovirus
<b>Target Names</b>	blaNDM-1
<b>Protein Names</b>	Recommended name: Beta-lactamase NDM-1 Short name= NDM-1 EC= 3.5.2.6Alternative name(s): Metallo-beta-lactamase NDM-1
<b>Expression Region</b>	29-270aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.