



Recombinant Rotavirus A Intermediate capsid protein VP6

Product Code	CSB-BP542165ROH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B1NKU0
Product Type	Recombinant Protein
Immunogen Species	Rotavirus A (strain RVA/Human/United Kingdom/ST3/1975/G4P2A[6]) (RV-A) (Rotavirus A (strain St. Thomas 3))
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MEVLYSLSKTLKDARDKIVEGTLYSNVSDLIQQFNQMIVTMNGNDFQTGGIGNL PIRNWTFDFGLLGTLLNLDANYVETARTTIEYFIDFIDNVCMDEMARESQRNG VAPQSEALRKLKLAGIKFKRINFNNSSEYIENWNLQNRRTGFVVFHKNIFPYS SFTLNRSQPMHDNLMGTMWLNAGSEIQVAGFDYSCALNAPANIQQFEHIVQL RRALTTATITLLPDAERFSFPRVINSADGATTWFFNPIILRPNNVEVEFLLNGQII NTYQARFGTIIARNFDTIRLSFQLMRPPNMPAVNALFPQAQPFQHHATVGLTL RIESAVCESVLADANETLLANVTAVRQEYAIPVGPVFPFGMNWTELITNYS EDNLQRVFTVASIRSMLIK
Source	Baculovirus
Protein Names	Recommended name: Intermediate capsid protein VP6
Expression Region	1-397aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.