



Recombinant Mouse Scavenger receptor cysteine-rich type 1 protein M130 (Cd163), partial

Product Code	CSB-EP650452MO-B
Abbreviation	Cd163
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q2VLH6
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	VVCQQLGCPTSICALGWANSSAGSGYIWMKVSCTGNESALWDCKHDGWG KHNTHEKDAGVTCSDGSNLEMRLVNSAGHRCLGRVEIKFQGKWGTVCDDN FSKDHASVICKQLGCGSAISFSGSAKLGAGSGPIWLDLACNGNESALWDCK HRGWGKHNCDAEDVGVICLEGADLSLRLVDGVSRCSGRLEVRFQGEWGT CDDNWDLRDASVVKQLGCPTAISAIGRVNASEGSGQIWLDNISCEGHEATL WECKHQEWGKHYCHHREDAGVTCS
Research Area	Others
Source	E.coli
Target Names	Cd163
Protein Names	Recommended name: Scavenger receptor cysteine-rich type 1 protein M130 Alternative name(s): CD_antigen= CD163 Cleaved into the following chain: 1. Soluble CD163 Short name= 2. sCD163
Expression Region	86-365aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.