



Recombinant Mouse Tyrosine-protein kinase JAK3 (Jak3), partial

Product Code	CSB-EP717252MO-B
Abbreviation	Jak3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q62137
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	LKYISLLGKGNFGSVELCRYDPLGDNTGPLVAVKQLQHSGPDQQRDFQREIQI LKALHSDFIVKYRGVSYGPGRQSLRLVMEYLPSCGLRDFLQRHRARLHTDRLL LFAWQICKGMEYLGARRCVHRDLAARNILVESEAHVKIADFGLAKLLPLGKDYY VVREPGQSPIFWYAPESLSDNIFSRQSDVWSFGVVLYELFTYCDKSCSPSAEF LRMMGPeregppLCRLLELLAEGRRLPPPPTCPTEVQELMQLCWAPSPHDP AFGTLSPQLDALWRGRPG
Research Area	Others
Source	E.coli
Target Names	Jak3
Protein Names	Recommended name: Tyrosine-protein kinase JAK3 EC= 2.7.10.2 Alternative name(s): Janus kinase 3 Short name= JAK-3
Expression Region	818-1100aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Target Details	This protein is a member of the Janus kinase (JAK) family of tyrosine kinases involved in cytokine receptor-mediated intracellular signal transduction. It is predominantly expressed in immune cells and transduces a signal in response to its activation via tyrosine phosphorylation by interleukin receptors. Mutations in this gene are associated with autosomal SCID (severe combined immunodeficiency disease).
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.