



Recombinant Human CMRF35-like molecule 1 (CD300LF), partial

Product Code	CSB-BP819470HU
Abbreviation	CD300LF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TDQ1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	TQITGPTTVNGLERGLSLTVQCVRSGWETYLKWWCRGAIWRDCKILVKTSGS EQEVKRDRVSIKDNQKNRTFTVTMEDLMKTDADTYWCGIEKTGNDLGVTVQV TIDPAPVTQEETSSSPTLTGHHLNDRHKLKLS
Source	Baculovirus
Target Names	CD300LF
Protein Names	Recommended name: CMRF35-like molecule 1 Short name= CLM-1Alternative name(s): CD300 antigen-like family member F Immune receptor expressed on myeloid cells 1 Short name= IREM-1 Immunoglobulin superfamily member 13 Short name=
Expression Region	20-156aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.