



Recombinant Human AP-2 complex subunit mu

Product Code	CSB-MP001872HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96CW1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MIGGLFIYNHKGEVLISRVRDDIGRNAVDVAFRVNVIHARQQVRSPVTNIARTSF FHVKRSNIWLAAVTKQNVNAAMVFEFLYKMCVMAAYFGKISEENIKNNFVLIY ELLDEILDFGYPQNSSETGALKTFITQQGIKSQHQTKEEQSQITSQVTGQIGWRR EGIKYRRNELFLDVLESVNLLMSPQQGQVLSAHVSGRVVMKSYLSGMPECKFG MNDKIVIEKQGKGTADETSKSGKQSI AIDDCTFHQCVRLSKFDSERSISFIPPDG EFELMRYRRTTKDIILPFRVIPLVREVGRTKLEVKVVIKSNFKPSLLAQKIEVRIPTP LNTSGVQVICMKGKAKYKASENAIVWKIKRMAGMKESQISAEIELLPTNDKKKK ARPPISMNFVFPFAPSGLKVRYLKVFEFPLNYS DHDVIKWVRYIGRSGIYETRC
Research Area	Transport
Source	Mammalian cell
Target Names	AP2M1
Expression Region	1-435aa
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length
Target Details	This gene encodes a subunit of the heterotetrameric coat assembly protein complex 2 (AP2), which belongs to the adaptor complexes medium subunits family. The encoded protein is required for the activity of a vacuolar ATPase, which is responsible for proton pumping occurring in the acidification of endosomes and lysosomes. The encoded protein may also play an important role in regulating the intracellular trafficking and function of CTLA-4 protein. Two transcript variants encoding different isoforms have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.