



Recombinant human Tubulin beta-2A chain

Product Code	CSB-EP025320HU-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q13885
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MREIVHIQAGQCQGNQIGAKFWEVISDEHGIDPTGSYHGDSDLQLERINVYYNE AAGNKYVPRAILVDLEPGTMDSVRSRGGPFGQIFRPDNFVFGQSGAGNNWAKGH YTEGAELVDSVLDVVRKESESCDCLQGFQLTHSLGGGTGSGMGTLISKIREE YPDRIMNTFSVMPSPKVSDTVVEPYNATLSVHQLVENTDETYSIDNEALYDICF RTLKLTPTTYGDLNHLVSATMSGVTTCLRFPGQLNADLRKLAVNMVFPRLHF FMPGFAPLTSRGSQQYRALTVPELTQQMFDSKNMMAACDPRHGRYLTVAAIF RGRMSMKEVDEQMLNVQKNSSYFVEWIPNNVKTAVCDIPPRGLKMSATFIG NSTAIQELFKRISEQFTAMFRRKAFHWTGEGMDEMEFTEAESNMNDLVSE YQQYQDATADEQGEFEEEEGEDEA
Research Area	Signal Transduction
Source	E.coli
Target Names	TUBB2A
Expression Region	1-445aa
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.