




Recombinant Rat Vitamin K-dependent protein S (Pros1)

Product Code	CSB-MP018754RA
Relevance	Anticoagulant plasma protein; it is a cofactor to activated protein C in the degradation of coagulation factors Va and VIIIa. It helps to prevent coagulation and stimulating fibrinolysis.
Abbreviation	Recombinant Rat Pros1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53813
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	ANTLLEETKKGNLERECIEELCNKEEAREVFENNPE TDYFYPKYL GCLGA FRV GAFSAARQSANAYPDLRSCVNAIPDQCDPMP C NEDGYLSCKDGGQGAFTCICK PGWQGD KCQFDINECKDPSNINGGCSQTC D NTPGSYHC SCKIGFAMLTNKKD CKDVDECSLKPSVCGTAVCKNIPGDFECECPNGYRYDPSSKSCKDVDECSEN TCAQLCVNYPGGYSCYCDGKKGFKLAQDQRSCEGIPVCLSLDL DKNYELLYLA EQFAGVVL YLKFRLPDITRFSAEFD FRTYDSEGI ILYAESLDHSNWLLIALREGKI EVQFKNEFSTQITTGGNVINNGIWNMVSVEELDDSVSIIKIAKEAVMNINKLGS LF KPTDGF LDTKIYFAGLPRKVESALIKPINPRLDGCIRGWNLMKQ GALGAKEIVE GKQNKHCFLTVEKGSYYPGSGIAQFSIDYNNVTNAEGWQINVT LNIRPSTGTG VMLALVSGDTPFALSLVDSGSGTSQDILVFVENSVA AHLEAITLCSEQPSQLK CNINRNGLELWTPVRKDV IYSKDLQRQLAILDKTMKGT VATYLG GVPDISFSAT PVNAFYSGCMEVNINGVQLDLDEAISKHNDIRAHSCPSVRKIQKNF
Research Area	Others
Source	Mammalian cell
Target Names	Pros1
Expression Region	42-675aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	73.9kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.