



Recombinant Human Talin-2 (TLN2), partial, partial

Product Code	CSB-MP897282HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y4G6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	RPQKIRMLDGSVKTVMVDDSKTVGELLVTICSRIGITNYEEYSLIQETIEEKKEE GTGTLKKDRTLRLDERKMEKCLKAKLHTDDDLNWL DHSRTFREQGV DENETLL LRRKFFYSDQNVDSRDPVQLNLLYVQARDDILNGSHPVSFEKACEFGGFQAQI QFGPHVEHKHKPGFLDLKEFLPKEYIKQRGAEKRFQEHKNCGEMSEIEAKVK YVKLARSLR TYGV SFFLVKEKMKGKNKLVPRLLGITKDSVMRVDEKTK EVLQE WPLTTVKRWAASPKSFTLDFGEYQESYYSVQTTEGEQISQLIAGYIDIILKKK
Research Area	Signal Transduction
Source	Mammalian cell
Target Names	TLN2
Expression Region	88-406aa
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Target Details	This gene encodes a protein related to talin 1, a cytoskeletal protein that plays a significant role in the assembly of actin filaments and in spreading and migration of various cell types, including fibroblasts and osteoclasts. This protein has a different pattern of expression compared to talin 1 but, like talin 1, is thought to associate with unique transmembrane receptors to form novel linkages between extracellular matrices and the actin cytoskeleton.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.