



Recombinant Human Multiple inositol polyphosphate phosphatase 1 (MINPP1)

Product Code	CSB-BP891977HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UNW1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	SLLEPRDPVASSLSPYFGTKTRYEDVNPVLLSGPEAPWRDPELLEGTCTPVQL VALIRHGTRYPTVKQIRKLRQLHGLLQARGSRDGGASSTGSRDLGAALADWPL WYADWMDGQLVEKGRQDMRQLALRLASLFPALFSRENYGRLRLITSSKHRC MDSSAAFLQGLWQHYPGLPPPVDVADMEFGPPTVNDKLMRFFDHCEKFLTE VEKNATALYHVEAFKTGPEMQNILKKVAATLQVPVNDLNADLIQVAFFTCSFDL AIKGVKSPWCDVFDIDDAKVLEYLNDLKQYWKRGYGYTINSRSSCTLFQDIFQ HLDKAVEQKQRSQPISSPVILQFGHAETLLPLLSLMGYFKDKEPLTAYNYKKQM HRKFRSGLIVPYASNLI FVLYHCENAKTPKEQFRVQMLLNEKVLPLAYSQETVS FYEDLKNHYKDILQSCQTSEECELARANSTSEDEL
Research Area	Signal Transduction
Source	Baculovirus
Target Names	MINPP1
Expression Region	31-487aa
Tag Info	Tag type will be determined during the manufacturing process.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.