



Rabbit anti-bovine IgG polyclonal Antibody, HRP conjugated

Product Code	CSB-PA00800F1Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	IgG purified from bovine serum
Raised In	rabbit
Species Reactivity	Bovine
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Alias	IgG
Product Type	Secondary Antibody
Immunogen Species	Bos taurus (Bovine)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.