



LIPG Antibody, Biotin conjugated

Product Code	CSB-PA897572LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5X9
Immunogen	Recombinant Human Endothelial lipase protein (360-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endothelial lipase (EC 3.1.1.3) (Endothelial cell-derived lipase) (EDL) (EL), LIPG
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	LIPG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.