



add Antibody

Product Code	CSB-PA001268LA01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22333
Immunogen	Recombinant Escherichia coli Adenosine deaminase protein (1-333AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adenosine deaminase (EC 3.5.4.4) (Adenosine aminohydrolase), add
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	add
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.