



ADAM12 Antibody, FITC conjugated

Product Code	CSB-PA001275LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43184
Immunogen	Recombinant Human Disintegrin and metalloproteinase domain-containing protein 12 protein (730-909AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Disintegrin and metalloproteinase domain-containing protein 12 (ADAM 12) (EC 3.4.24.-) (Meltrin-alpha), ADAM12, MLTN
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ADAM12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.