



AKR1B1 Antibody, Biotin conjugated

Product Code	CSB-PA001539LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15121
Immunogen	Recombinant Human Aldose reductase protein (114-230AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aldose reductase (AR) (EC 1.1.1.21) (Aldehyde reductase) (Aldo-keto reductase family 1 member B1), AKR1B1, ALDR1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AKR1B1