




ASS1 Antibody

Product Code	CSB-PA002234LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00966
Immunogen	Recombinant Human Argininosuccinate synthase protein (168-222AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:2000-1:10000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Argininosuccinate synthase (EC 6.3.4.5) (Citrulline--aspartate ligase), ASS1, ASS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ASS1

Image



Western Blot

Positive WB detected in: Mouse brain tissue,

Hela whole cell lysate


All lanes: ASS1 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 47 kDa

Observed band size: 47 kDa



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA002234LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.