



ASS1 Antibody, HRP conjugated

Product Code	CSB-PA002234LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00966
Immunogen	Recombinant Human Argininosuccinate synthase protein (168-222AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Argininosuccinate synthase (EC 6.3.4.5) (Citrulline--aspartate ligase), ASS1, ASS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ASS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.